CENTRAL VENOUS ACCESS CATHETERS

Chad Bird, B.S. NREMT-P

Objectives

- What are Central Venous Access Catheters (CVAC)
- Who will have a Central Venous Access Catheter
- Types of Central Venous Access Catheters
- Accessing the Central Venous Access Catheters
- Complications with Central Venous Access Catheters

What is a CVAC?



http://www.med.unsw.edu.au/SPHCMWeb.nsf/pageprintfriendly/ResStrengthInfectious?opendocument

What is a CVAC?

- A tube is inserted directly into a central vein with the distal tip placed just above the right atria.
- Common veins include but are not limited to:
 - External Jugular
 - Subclavian
 - Very popular due to decreased risk of infection
 - Brachial
 - Femoral
 - Has a higher risk of infection




Figure 42-20 A, Central venous catheter insertion and exit site. (Aehlert, Barbara. *Paramedic Practice Today: Above and Beyond.* Mosby, 022009.).

Who might have a CVAC?

- A CVAC can minimize painful peripheral IV attempts, so patients requiring frequent IV access may be candidates.
- Other reasons for insertion of a CVAC are:
 - Patient has poor peripheral venous access.
 - Patient requires medications that are given routinely or that irritate peripheral veins.
 - Patient requires rapid fluid or medication administration.

Patients to suspect of having a CVAC include, but are not limited to ...

Patients to suspect of having a CVAC include, but are not limited to ...

Cancer patients

Patients to suspect of having a CVAC include, but are not limited to ...

- Cancer patients
- Special needs patients

Patients to suspect of having a CVAC include, but are not limited to ...

- Cancer patients
- Special needs patients
- Dialysis patients

Patients to suspect of having a CVAC include, but are not limited to ...

- Cancer patients
- Special needs patients
- Dialysis patients
- Critical care patients being transferred from one facility to another

Types of CVAC


- Tunneled catheters are the most common. They are tunneled under the skin into a central vein.
 - Hickman®, Broviac®, or Groshong® catheters
 - May have one, two, or most commonly three lumens
 - Lumens may be used for fluids and meds, drawing blood, invasive monitoring, or nutritional support
 - These will typically exit the chest making them easy to miss if a physical inspection is not performed

Types of CVAC A C GROSHONG DOUBLE (II ATMILE) C FIGURE 21-3 A, Hickman catheter. B, Broviac catheter. C, Groshong catheter.

(From Clayton BD, Stock YN: Basic pharmacology for nurses, ed 13, St Louis, 2004, Mosby. Courtesy Chuck Dresner.)
(Ogden, Sheila J. Calculation of Drug Dosages, 8th Edition. C.V. Mosby, 032007.).

Types of CVAC, cont.

■ PICC (peripheral inserted central catheter) lines are inserted though a peripheral vein such as the antecubital fossa rather than the chest or neck, but the distal tip is still placed just above the right atria.



http://terisblip.com/2010/10/18/mn-picc-thanks-force-blog-lag-good-night/

Types of CVAC, cont.

- Implanted ports, also know as a port-a-cath, are similar to tunneled catheters but do not exit the skin.
 - Small reservoir implanted under the skin in the upper anterior chest
 - Greatly reduce the risk of complication
 - Can be very easy to miss on physical inspection.

Types of CVAC, cont.



FIGURE 21-5 Silicone venous catheters with infusion ports. (From Clayton BD, Stock YN: Basic pharmacology for nurses, ed 13, St Louis, 2004, Mosby. Courtesy Chuck Dresner.)

(Ogden, Sheila J.. Calculation of Drug Dosages, 8th Edition. C.V. Mosby, 032007.)

Accessing a CVAC

- Gather the equipment.
 - 12 cc syringe or preferably a 20 cc syringe.
 - Multiple syringes if you are going to draw blood.
 - Providone iodine or alcohol swabs.
 - Huber needles for implanted devices.
- The most important thing to remember is to use aseptic techniques. Central lines are susceptible to infections which can be fatal.
 - Ensure proper BSI precautions.
 - $\ ^{\square}$ Gloves, eye protection, and mask
 - If possible, wash hands prior to accessing the CVAC.
 - Providone iodine should be used when possible.

Accessing a PICC line

- To access a PICC catheter:
 - BSI!!!
 - Swap with providone iodine or alcohol wipes.
 - Without touching the end of the catheter, attach syringe to cap or insert the needle.
 - Unlock the clamp.
 - Gently aspirate at least 10 mL of blood.
 - Gently flush the line with 10 mL of NaCl.
 If line does not flush with ease, do not use.
 - If line flushed without problem, you may begin the infusion.

Accessing a Tunneled Catheter


- Accessing a Hickman®, or Broviac® catheter is essentially the same procedure as accessing a PICC line with a couple of points to remember:
 - Typically these catheters are flushed with heparin, so be sure to aspirate at least 10 mL of blood.
 - Remember to close the clamps when not aspirating or flushing the catheter due to lack of anti-reflux valves.


Accessing an implanted device

- The process for accessing an implanted device will differ slightly.
 - BSI!!!
 - If possible, wash your hands and utilize sterile gloves.
 - Do not delay access in critical situations. If you can not use sterile procedures, utilize clean gloves.
 - Swab the skin with providone iodine or alcohol.
 - Access the port using a Huber needle. Stabilize the port with your thumb and index finger. Inserting the needle at a 90 degree to the skin until it stops.
 - You may need to pull the needles back slightly if unable to aspirate blood.

Accessing an implanted device, cont.

- Once inserted, place sterile gauze under the wings of the needle and secure using an occlusive dressing or tape.
- Since these ports are often heparinized, it is extremely important to aspirate 10 mL of blood prior to flushing or administering a bolus.
- Obtain blood samples if needed and your ready to flush and administer fluids and/or medication.
 - Be sure to clamp any tubing when not in use.





Complications with CVAC

- □ Central catheters are relatively safe for paramedics to access in the clinical setting or in pre-hospital setting, but complications may occur. Proper training and remaining calm will help you troubleshoot the problem and determine a solution.
- Prevent the complication.
- Some complications may be fatal!!!

Complications with CVAC, cont.

- Infections
 - Infections are not always preventable, but using sterile or clean techniques can help avoid contamination.
 - Wash hands if possible.
 - Always use BSI precautions.
 - Clean the port using aseptic techniques any time you access for blood draws, meds, or fluid administration

Complications with CVAC, cont.



http://blog.timesunion.com/mdtobe/category/legal-issues/

Complications with CVAC, cont.

- Poor flow
 - The CVAC may have poor flow. One possible cause could be that the line is positional. Try having the patient turn their head to one side and/or take a deep breath.
- Poor flow could indicate a more serious complication such as microthrombi.
 - Important to aspirate 10 mL prior to flushing or infusing meds to prevent flushing micro clots into the vasculature.

Complications with CVAC, cont.

- Bleeding around tube can be common and is typically not life threatening if addressed early.
 - This can indicate catheter dislodgment.
 - Try to avoid dislodgment by not letting the tube get snagged on bed rails or other objects.
 - Treatment in the field is bleeding control.

Summary

- EMS provider may encounter one of many types of central venous access catheters in a clinical or pre-hospital setting.
- Familiarize your self with these devices and the patients who may have them, and be prepared to access them in an emergency situation.
- Proper training and practice will increase the chance of proper use and minimize the chance of potentially life threatening complications.
 - Google, YouTube, ect.

References

- □ Pollak, Md, Faaos, A. N. (Ed.). (n.d.). *Critical Care Transport. Sudbury, MA: Jones & Bartlett.*
- ☐ John's Hopkins University (2011). www.hopkinsmedicine.org
- Vanderbilt University Medical Center (2011). www.mc.vanderbilt.edu
- Ems World (January 12, 2011).